HYA CORP *endo gel*

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Datum der Information: 05.2010

Product information

Post surgical tissue adhesion is a common consequence of surgical procedures, these post-surgical adhesions are due to the formation of areas of contact, made of fibrous tissue, between adjacent internal organs, which may induce different post operative symptoms like: pain, infertility, bowel obstruction, restriction of movement in case of tendon and spinal surgery, and thus a second surgery for removing fibrous tissue would be needed.

In accordance to excellent surgery practice and preventing formation of post surgical adhesions it is recommended to use

HYAcorp endo gel; which forms a physical barrier against the contact between adjacent tissues and remain on the site of application for seven to ten days, which is sufficient to avoid the formation of adhesions during the critical healing process of injured tissues.

HYAcorp endo gel is a bio-degradable, sterile, transparent, high viscous gel obtained by crosslinking special designed hyaluronic acid, one of the main components of human connective tissue and of epithelial and mesothelial tissues and it perfectly adheres to the tissue surface and to the abdominal wall creating an anti adhesion barrier.

Indications

HYAcorp endo gel is an absorbable, anti-Adhesion barrier gel that helps to reduce the post surgical adhesion formation in patients undergoing laparotomic or laparoscopic abdominal, pelvic-gynecological and orthopedic surgeries.

Instructions for use

- 1 - After having brought the product to room temperature, open the pouch and transfer the syringe into the sterile field, adopting the normal aseptic techniques used in the surgical theatres.
- 2 - Remove the protective cap on the tip of the syringe and insert the special cannula/applicator.
- 3 - After aspiration of all irrigation fluid, apply the gel inside to completely cover all anatomical sites where adhesion prevention is desired, even if the thickness of the gel applied does not influence the efficacy of the product, it is recommended to apply a 1-2 mm thick layer of gel.
- 4 - Do not irrigate after application of the product.

Special instructions for pelvic-gynecological surgeries

- 1 - Lift ovary away from pelvic sidewall and apply a single layer of gel to cover the ovarian fossa and posterior surface of the ovary.
- 2 - Return ovary to normal anatomical position and apply a single layer of gel to cover the anterior portion of the ovary.
- 3 - Apply a single layer of gel to cover the fallopian tube, including the ampulla and mesosalpinx.
- 4 - Apply a single layer of gel to cover the lateral aspect of the uterus facing the adnexa.

After placement - all surgeries:

- 1 - Do not apply irrigation fluid to the sites coated with **HYAcorp endo gel**.
- 2 - Do not attempt to reposition gel once applied.
- 3 - After completely covering the site once, do not apply additional gel to that site, if the site is not adequately covered by a single layer of gel, additional layers can be applied next to previously applied gel. Avoid applying additional layers of gel on top of one another.

Contraindications

Known hypersensitivity to hyaluronic acid and Gram positive bacteria. As all other implant materials; the device must not be used in patients with active infection and/or contamination of the surgical site.

Storage

Storage temperature 2 °C – 25 °C.
Do not freeze.

Warning and precautions

This product is to be used only under surgical aseptic conditions by a licensed medical practitioner. Once package is opened discard any remaining gel after use. Do not use if package is damaged.

Presentation

The product is designed to satisfy all needs of different surgical procedures utilized by Surgeons, Gynecologists and Orthopedics:

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D Gebrauchsanweisung beachten
GB Refer to instruction for use

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BioScience GmbH