# EPHYAL®

D

## GEBRAUCHSINFORMATION

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### WARNHINWEIS

### ANWENDUNG

### LIEFERFORM

### HERSTELLER

### SYMBOLE

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## USE INSTRUCTIONS

### COMPOSITION

**EPHYAL** is a viscoelastic, isotonic, sterile solution. It contains hyaluronic acid sodium salt of non animal origin. Hyaluronic acid has a low molecular weight (VLMW) which guarantees excellent absorption in the epidermis and the dermis.

### AREAS OF APPLICATION

Hydration of tissue. To be applied to the forehead, to the area around the eyes and cheeks as well as to the area of the nose and mouth, to the lips, the neck, the lower neckline and to the hands.

### CONTRAINDICATIONS

Do not use **EPHYAL** if you are known to suffer from hypersensitivity towards hyaluronic acids.

### INTERACTIONS

**EPHYAL** is not so far known to be incompatible with other topically applied products.

### FURTHER INFORMATION

**EPHYAL** is only suitable for external use. Do not use damaged ampoules. Do not use product after the shown expiry date. Keep product in a place which is not accessible to children.

### APPLICATION

In order to hydrate the face, neck and the lower neck areas, apply one ampoule of **EPHYAL** once a day on previously washed skin. If required (for extremely dry and damaged skin) **EPHYAL** can also be applied several times a day. **EPHYAL** penetrates the skin and achieves a silky effect very quickly.

In contrast to cosmetics, which according to the currently regulations, are only able to work on the surface of the skin, **EPHYAL** is an authorised, sterile medicinal product. It permeates the lower layers of the skin including the dermis, which is primarily responsible for the formation of wrinkles.

### FORM OF SUPPLY

50 ml vial containing 40 mg hyaluronic acid (VLMW)

### MANUFACTURER

BioScience GmbH  
Friedenstr. 29  
56427 Siershahn  
Germany

### SYMBOLS

Sterilisation with steam

Not for reuse

Pay attention to instructions for use

Storage temperature

Usable until

Batch designation

Date of manufacture

CE0297

CE number commensurate with MDD 93/42/EEC  
0297 is the number of the „denominated location“

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